

Experiment 4: Gravimetric Determination of Lead Chromate

SYNOPSIS Lead is measured by precipitation with chromate, dried, and weighed.

READINGS: Pages 269-273 in Critical Reviews

SOLUTIONS

Concentrated NaOH

pH meter

0.1% Nitric Acid (Conc. $\text{HNO}_3 = 76\%$, so about 1 mL conc./100 mL water)

0.10 M chromium nitrate, $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ (4 g/100 mL)

0.12 M potassium bromate, KBrO_3 (2 g/100 mL)

Acetate buffer solution: 6 M in acetic acid, 0.6 M in sodium acetate.

GLASSWARE

2 erlenmeyer flasks

2 porous porcelain filter crucibles

suction adapters

PROCEDURE

1. To each flask add 4 mL of 1000 ppm Pb. Bring to 20 mL vol.. To the fifth flask add 20 mL of one of your soil sample digests. To one of the first four (lead standard) flasks add 4 mL of 1000 ppm Zn.
2. If necessary, neutralize 20 mL solution with NaOH to pH 7 (use indicator paper) in a 100 mL beaker. Solution will be slightly cloudy. *What is the chemistry that makes it cloudy?*
3. Adjust volume of sample to about 20 mL.
4. Add 10 mL $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ solution and 10 mL of KBrO_3 . Solution will be clear blue.
5. Heat but do not boil for 30 minutes.
6. When solution is clear and yellow (= measure of extent of chromic oxidation), add 10 mL of buffer and heat 5 more minutes.
7. Weigh your 6 filter crucibles.
8. Cool the mixture and filter off the lead chromate on a sintered glass or porous porcelain filter crucible.
9. Wash the precipitate with 2 or 3 small portions of 0.1% nitric acid.
10. Dry at 120°C for 30 minutes.
11. Cool and weigh as PbCrO_4 ($\text{Pb}/\text{PbCrO}_4 = 0.641108$).

REPORT

1. What is the chemistry of all of the reactions involved in forming lead dichromate? Write out

balanced reactions.

2. What is the chemistry that makes it cloudy?
3. What is the relationship between the solubility of PbCrO_4 and its use as a yellow pigment for traffic paints?
4. What is the purpose of the added KBrO_3 ?
5. What is the relative standard deviation of your 3 ppm sample?
6. Does the addition of Zn to the standard affect your analysis?
7. What is the absolute error of your 3 ppm sample? What might cause the error?
8. How easy would it be to explain this method to a technician or a 5th grader?