



~~_____~~
~~_____~~
~~_____~~

Chemistry 311
November 22, 1996

Arrangement of Responsibility in Report Preparation for Experiment 25: ICP-MS

Data Collection

~~_____~~
~~_____~~

Analysis and Manipulation of Data

~~_____~~
~~_____~~
~~_____~~

Majority of Report Writing

~~_____~~

Editing

~~_____~~
~~_____~~

Signed:

~~_____~~


~~_____~~

~~_____~~

~~_____~~

~~_____~~

~~_____~~


November 22, 1996

Experiment 25: ICP-MS Determination of Lead

I. Introduction

The Inductively Coupled Plasma - Mass Spectroscopy was chosen as the method for determining the amount of lead in the samples of house dust obtained by the schoolchildren. This method was decided upon for many reasons which include a low limit of detection, an instrument that is easy to use, and the little amount of time needed to analyze many samples. Because the instrument is a mass spectrometer, the system can be adjusted to measure the amounts of selected species such as the various isotopes of lead (^{204}Pb , ^{206}Pb , ^{207}Pb , and ^{208}Pb) along with internal standards (^{50}V , ^{51}V , ^{121}Sb , ^{123}Sb , and ^{232}Th).

II. Materials and Methods

Solutions: 70 % HNO_3 stock solution
 1000 ppm Vanadium stock solution
 1000 ppm Antimony stock solution
 1000 ppm Thorium stock solution

Prepared Solutions: 1 % HNO_3 - wash solution
 0.1 % HNO_3 with 50 ppb each of V, Sb, Th - Internal standards and blank
 Pb standard solutions (10, 50, 100, 200, 500, and 1000 ppb)
 Sample solutions of internal controls (diluted x 5)
 Samples of digested baby wipes (diluted x 5)

Materials: 18 50 ml volumetric flasks
 2 1 L volumetric flasks
 micropipettes
 10 ml graduated cylinder
 Parafilm

Instrumentation: Inductively Coupled Plasma - Mass Spectrometer (ICP-MS)

Setup:	Element	Scanned Regions
	V	50, 51
	Sb	121, 123
	Pb	204, 206, 207, 208
	Th	232

Procedure: Pumping → Standby → Operating Mode

Methods: For preparing the solutions, Sabeen and Kim made up the wash solution and the solution containing the internal standards. The wash solution was made by adding 14.3 ml of the 70% HNO₃ stock solution to a 1 L volumetric flask that contained a little bit of deionized water. The solution was then diluted to 1 L with deionized water. The 0.1% HNO₃ solution was made with 1.43 ml of the HNO₃ stock solution and 50 ml each of the V, Sb, and Th stock solutions which were diluted to 1 L in a 1 L volumetric flask with deionized water. June made up the lead standard solutions by micropipetting 5, 25, 50, 100, 250, and 500 µl of 100 ppm Pb stock solution into six 50 ml volumetric flasks. Sabeen diluted all of the Pb solutions with the 0.1% HNO₃ solution while Kim covered them with parafilm and shook the flasks. June also measured out 10 ml of the lead solutions from the digestion and poured them into the remaining 50 ml volumetric flasks. Sabeen diluted them with the 0.1% HNO₃ solution, and Kim made sure that each solution was thoroughly mixed by shaking the flasks. In analyzing the solutions, Kim ran the computer while June transferred the intake tube from the wash solution to the various samples.

III. Results

The raw data obtained are contained in the appendix. Otherwise, the data are summarized in the following table (Table 1) and in Table 2 and Figure 1 which are on one of the following pages.

Table 1: The amount of lead in the samples of KimWipes and house dust using the data on ²⁰⁶Pb.

Sample	Concentration (ppm)	Concentration (µg/ft ²)
Sabeen's internal control 1	4.350	435.0
Sabeen's internal control 2	4.530	453.0
June's internal control 1	4.710	471.0
June's internal control 2	5.025	502.5
Kim's internal control 1	4.440	444.0
Kim's internal control 2	4.470	447.0
Sabeen's Unwashed Sample	1.360	136.0
Sabeen's Washed Sample	0.1975	19.75
June's Unwashed Sample	0.4140	41.40
June's Washed Sample	0.5600	56.00
Kim's Unwashed Sample	0.0125	1.250
Kim's Washed Sample	0.0000	0.000

The internal controls are the digestions performed with 5.0 ppm of lead added to a KimWipe in a beaker. These digestions were performed twice by each student so they were given arbitrary numbers of 1 or 2 for each student.

Kim's Washed sample had counts with negative values so the concentration was assumed to be 0.0 ppm. Her Unwashed sample had negative values for ²⁰⁴Pb and ²⁰⁷Pb so the concentration of lead is probably near zero.

The concentrations already took the 5 times dilution into account which means that the raw data was multiplied by 5 to obtain the concentrations in ppm. The concentrations in ppm were then used to calculate the amount of lead in µg/ft².

IV. Summary/Discussion

According to Table 1, the children's samples are definitely below the HUD limit of $500 \mu\text{g}/\text{ft}^2$. However, some error could be attributed to the children who might not have performed their part of the collection of the samples well. They might not have measured a square foot perfectly or understood the instructions concerning collecting the sample of an unwashed sample versus a washed sample particularly since June's washed sample had a higher lead content than the unwashed surface and was a noticeably dirtier baby wipe than the baby wipe for the unwashed surface before the digestions. If Sabeen's samples were acquired properly, they proved that lead content could decrease significantly with washing.

Another cause of error would be the digestions themselves which were difficult to perform. The error can be shown by the measurements of the internal standards. These internal standards originally had 5.0 ppm lead, but were analyzed with the ICP-MS to have less than this amount for everyone's samples except for June's second one which exceeded the original concentration. This may have also been caused by the digestion process which took a while and was altered. Another indication of human error is the counts per second of vanadium, antimony, and thorium in Table 2 for the controls and the samples. Each of those solutions were to contain the same amount of internal standards because they were the 5 times dilutions of the samples and thus had 10 ml each of the sample and 40 ml each of the 0.1 % HNO_3 with the internal standards. However, the counts per second of each of the elements varies significantly among the samples even though they should be the same. One reason for the inconsistency is that the 10 ml of samples were not measured accurately or the dilutions with deionized water was not done accurately. This shows the great probability of human error.

The instrument itself could also have been a cause of error which can be exemplified by the limit of detection (LOD) and the limit of quantification (LOQ) but also includes human error in the ability to make the lead standard solutions. Since none of the solutions were run repeatedly, the standard deviation of the blank would have to be taken from the regression line as being the same as the mean of the blank. From Figure 1, the mean of the blank as well as the standard deviation of the blank was determined to be 461.547 counts per second. Since LOD equals the mean of the blank plus three times its standard deviation, the LOD was calculated to be 1846.188 counts per second ($461.547 + 3 \times 461.547$) which pertains to 14.35 ppb Pb. The LOQ is the mean of the blank plus nine times its standard deviation which would be 4615.47 counts per second. This pertains to 43.05 ppb Pb which is a very small number and means that the instrument can make very accurate measurements. Furthermore, the linear range reached 1000 ppb Pb but could extend further since the 1000 ppb was the maximum concentration that we used. Finally, the r value was determined to be 0.99982 which indicates that the regression line is a very good fit.

The relatively little background is one reason for the low LOD, and the ICP-MS has many causes for such little background noise. One of these is the use of argon plasma as an ionizing source. At the temperature of the plasma, Ar is only slightly ionized (0.002 %) so any ions of Ar will not contribute significantly to the detector and cause much background. Additionally, the mass of Ar (40) is far from any of the isotopes being analyzed so it will not interfere. Another contributor to the little background is the high vacuum behind the skimmer. The decrease in pressure increases the mean free path which means that the ions travel separately instead of colliding with one another. This allows the ions to hit the detector one at a time so they will not interfere with each other.

The primary reason why the ICP-MS has such a low LOD is the efficiency of the instrument. The instrument can make many measurements in a short amount of time. In 50 μ s, about 10^7 counts can be acquired. With a coadded spectra, the signal-to-noise ratio would increase so much that any noise would be practically insignificant.

The Gaussian shape curve could be attributed to the different energies of the ions as they travel through the quadrapole. The ions cannot enter the quadrapole with the exact same velocity. In fact, the velocity at which the ions enter are random which is shown as the Gaussian curve. The curve does not mean that the ions have varying masses because they do not. In fact, each ion of the same isotope should have the same mass because they have the same composition of neutrons and protons that determine contributes to its mass. Perhaps, this would be the main instrumental error because the peaks are not perfectly narrow. Every other part of the instrument would give narrow peaks which would have a low degree of uncertainty.

Since the first ionization energy of lead is 7.4 eV and the percent ionized at a plasma temperature of 7500 K is 0.994 %, the formation of Pb^+ would not be a problem. In fact, 97.5 of Pb is completely ionized which would mean that the singly ionized lead would not interfere significantly with the measurements.

The excellent choice for using this method to analyze the samples is supported by many reasons. First of all is the low LOD and LOQ along with the large linear range and the good r value. Another reason for making ICP-MS the primary method is the short turn around time for measuring samples which was approximately 5 minutes. This method is also ideal because of the lack of use of hazardous materials except for the lead that is analyzed. Moreover, a paper trail could be easily constructed for this procedure because it does not require much preparation. Unfortunately, one flaw is the difficulty in instructing a technician on how to use the instrument. In fact, delays were caused because a few times someone who understood the computer program were not available so the experiment could be run. However, with a knowledgeable technician nearby, the procedure was simple and accurate.

Instrument	LOD (ppb)	LOQ (ppb)
UV-Vis	2400.0	indeterminable
GFAA	169.0	508.0
Fluorescence	2290.0	2360.0
ASV	31.3	93.9
ICP-MS	14.35	43.05

V. Sample Calculations

Determination of $\mu\text{g}/\text{ft}^2$:

Sabeen's Unwashed Sample

$$1.360 \text{ ppm} = \frac{1.360 \mu\text{g Pb}}{1 \text{ ml digestion}} \times \frac{100 \text{ ml digestion}}{1 \text{ ft}^2} = 136.0 \mu\text{g}/\text{ft}^2$$

ICP-MS

Start Mass : 46.6
 End Mass : 234.4
 Mass Calibration : $y = a + bx + cx^2 + \dots + ex^4$
 -1.37E-01 4.44E-03 -4.57E-10 0.00E+00 0.00E+00

TABLE 2

Blank - 1% Nitric Acid

Mass	% Ab	Z	Int CPS
49.95	0.24	23	37.71
50.94	99.76	23	58.73
120.9	57.25	51	176.19
122.9	42.75	51	127.97
203.97	1.37	82	65.53
205.97	25.15	82	1726.77
206.98	21.11	82	2575.6
207.98	52.38	82	5299.47
232.04	100	90	34.62

**Internal Standard - 0.1% Nitric Acid
with 50 ppb V, Sb, Th**

Mass	% Ab	Z	Int CPS
49.95	0.24	23	207.72
50.94	99.76	23	42010.35
120.9	57.25	51	15818.22
122.9	42.75	51	10547.15
203.97	1.37	82	7.42
205.97	25.15	82	24.73
206.98	21.11	82	25.96
207.98	52.38	82	55.64
232.04	100	90	52220.88

10 ppb Pb

Mass	% Ab	Z	Int CPS
49.95	0.24	23	192.26
50.94	99.76	23	40909.35
120.9	57.25	51	14946.31
122.9	42.75	51	10124.63
203.97	1.37	82	73.57
205.97	25.15	82	1354.57
206.98	21.11	82	2432.78
207.98	52.38	82	4698.43
232.04	100	90	50884.63

50 ppb Pb

Mass	% Ab	Z	Int CPS
49.95	0.24	23	152.7
50.94	99.76	23	39743.03
120.9	57.25	51	14448.8
122.9	42.75	51	9851.19
203.97	1.37	82	301.07
205.97	25.15	82	7379.8
206.98	21.11	82	12850.12
207.98	52.38	82	25146.49
232.04	100	90	48363.19

Multi-Element Calibration Fri Nov 22 1996

Calibration graph coefficients for selected isotopes.

User Name : STUDENT Experiment name : LEAD2

Analysis procedure : LEAD STANDARD and samples

Introduction method : 60 sec. acquire method 30up 10wash

First sample started at : Fri Nov 22 1996 11:31:03

Last sample started at : Fri Nov 22 1996 12:14:23


Polynomial fit : $y = a_0 + a_1 \cdot x + a_2 \cdot x^2$

Units of response : $y =$ counts per second

Units of concentration : $x =$ ppb

<u>Element</u>	<u>Symbol</u>	<u>Mass</u>	<u>a0</u>	<u>a1</u>	<u>a2</u>	<u>regression</u>
Lead	Pb	204	51.6	3.80	0.0	1.000
Lead	Pb	206	462	96.5	0.0	1.000
Lead	Pb	207	1128	166	0.0	1.000
Lead	Pb	208	1614	330	0.0	1.000


Pb 204



Coefficients : $Y = A_0 + A_1X + A_2X^2 \dots A_nX^n$

A0	A1	A2	A3	A4	Regression
5.15608E+1	3.79688E+0	0.00000E+0	0.00000E+0	0.00000E+0	9.99666E-1

sabeen washed 1



EXPERIMENT 05: ICP-MIC

Nov 20, 2016

RUN AT 5100

Pb STANDARDS in 50ml volumetric

Pb = 100ppm	10 PPB	50ul
	50	250ul
	100	500ul
	200	1000ul
	500	2500ul
	1000	5000ul

INTERNAL STANDARDS in 1L volumetric

.1% HNO₃ (70%)

$$V_1 = \frac{(1000)(10000\text{ml})}{100000} = 100\text{ml}$$

.1
(100ppm)

50 PPB V, Sb, Th = 50ul of each

RUNNING SAMPLES

1. Blank 1% HNO₃
2. 1% HNO₃ - V, Th, Sb - INTERNAL STANDARD
3. 10 PPB Pb
4. 50 PPB Pb
5. 100 PPB Pb
6. 200 PPB Pb
7. 500 PPB Pb
8. 1000 PPB Pb
9. Kim control 1
10. Sabreen control 1
11. June control 1
12. Kim control 2
13. Sabreen control 2
14. June control 2
15. Sabreen unwashed
16. Sabreen washed
17. Kim unwashed
18. Kim washed
19. June unwashed
20. June washed

4/24/10

ICP-MS - used a known amount of Pb in hows dust
stock solutions: 70% HNO₃
1000 ppm V, Sb, Th each
100 ppm Pb

Prepared solutions:

0.1% HNO₃ with 50 ppb each of V, Sb, Th - Internal Standards

To a 1L volumetric flask containing some DI water, 1.43 ml of the 70% HNO₃ stock solution was added with a volumetric pipette by Sigma. 50 µl each of V, Sb, and Th were subsequently micro-pipetted. The solution was diluted and shaken by hand.

$$M_1 V_1 = M_2 V_2$$

$$V_1 = \frac{M_2 V_2}{M_1} = \frac{(0.1\%)(1000 \text{ ml})}{70\%} = 1.43 \text{ ml}$$

1.0% HNO₃ waste solution

143 ml of 70% HNO₃ was diluted to 1L in a 1L volumetric flask.

All of the glassware was rinsed with diluted nitric acid, then, and deionized water prior to use.

11/22/96

CP-713

samples + samples : 5 - dilution (10 ml in 50 ml vol flasks)

Conc standards from 100 ppm Pb stock sol'n

Conc (ppm)	Vol of stock sol'n (ml)	to 50 ml vol flasks
10	5	
50	25	
100	50	
200	100	
500	250	
1000	500	

Kim measured. Seven dilutes w/ Int. Std., Kim covered w/ parafilm + shook the

Sequence of analysis

1. 0.1% conc. V, Se, Th
2. 10 ppm Pb
3. 50 ppm Pb
4. 100 ppm Pb
5. 200 ppm Pb
6. 500 ppm Pb
7. 1000 ppm Pb
8. Kim control 1
9. Seven control 2
10. Kim control 1
11. Kim control 2
12. Seven control 2
13. June control 2
14. Seven control 2
15. Seven control 2
16. Kim blank
17. Kim blank
18. Kim blank
19. Kim blank

The dilutions, June measured 10 ml with a graduated cylinder. Seven dilutes with Internal Std, Kim covered with parafilm and shook.

~~XXXXXXXXXX~~
ICP-MS Second Run

11/22/11

dilution of samples five fold
ex. 10ml of sample solution
brought up to 50ml with
addition of 0.1% HNO₃ w/ Internal
Standards

Data Lead Samples

1. blank 0.1% HNO₃
2. 0.1% HNO₃ plus Internal Standards
3. 10ppb Pb
4. 50ppb Pb
5. 100ppb Pb
6. 200ppb Pb
7. 500ppb Pb
8. 1000ppb Pb
9. Kimcontrol 1
10. June control 1
11. Sabien Control 1
12. Kimcontrol 2
13. June Control 2
14. Sabien Control 2
15. Sabien unwashed
16. Sabien washed
17. Kim unwashed
18. Kim washed
19. June unwashed
20. June washed